

Biochemistry/
selection
↓

① Shotgun Sequencing

WGS
RNA seq
shotgun metagenome
ChIP seq

Sample
(DNA/RNA)

→ (All DNA → Fragmented)

Digitized

↓
Computer

1. Random sampling
2. Quantitative


Not Shotgun:
IGS / barcode
tag end
exome
RAD seq

Illumina


DNA sample
(or cDNA)

→ fragment ~ 300-800 bp
"single molecule"

== == ==
ligate
primers



PCR
BAC
bias against GC/AT
rich




x 20m - 2 billion

sequence
in situ

2m - 2bn
sequences

→ 100-300 bp each

What do we want? (Interlude)



Concerns:

- ① we want haplotype resolution - which variants are on the same physical chr.
- ② we want high coverage - we want to see all PstI at least 5-6 times

latency

throughput
(bp per time)

cost

accessibility /
realtime

investment cost


Long read sequencing

- \$ but a lot longer than Illumina
(~ 2 kb - 100 kb)
- single molecule (no PCR)
 - more input PstI
 - lower throughput

Nanopore + Pacific Biosciences

Two reasons:

- ① single molecule → PCR biases, don't use !!
 - ② long read
 - ↳ haplotype resolution
 - ↳ repeat resolution
 - ↳ good for structural variants
 - long indels
 - inversions
 - etc.
- no amplify
↳ minimal bias



Other:

- RNAseq
- Metagenomics

